ENG COMPRESSOR NEBULIZER LD. *Instruction Manual.*
ROU INHALATOR COMPRESOROWY LD. *Instrucția de utilizare.*
HUN COMPRESSZOROS INHALÁTOR LD. *Szórószagt. Használati útmutató.*
ROU APARAT DE AEROSOLI CU COMPRESOR LD. *Manual de instrucțiuni.*
BGH КОМПРЕСОРЕН ИНХАЛАТОР LD. *Ръководство за експлоатация.*

fig.1 rys.1 1. abra fig.1 рис.1
 PARTS AND COMPONENTS DENUMIREA PĂRȚILOR ȘI COMPONENTELOR
 PODSTAWOWE CZĘŚCI I KOMPONENTY НАЗВАНИЯ НА ЧАСТИ И КОМПОНЕНТИ
 ALKATRESZEK MEGNEVEZÉSE

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
fig.3 rys.3 3. abra fig.3 рис.3
 Depending on the type of baffle, the particles of different sizes are distributed aerosols as follows:

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ENG

- PARTS AND COMPONENTS (fig.1)**
- COMPRESSOR.
 - POWER SWITCH.
 - ANGULAR HOLDER.
 - CONNECTOR OF COMPRESSOR.
 - INTAKE.
 - SOCKET FOR THE FILTER.
 - NEBULIZER LD-N105.
 - UPPER PART OF NEBULIZER.
 - BAFFLE "B" (BLUE) LD-N002.
 - BAFFLE "C" (RED) LD-N002.
 - NOZZLE.
 - BOTTOM OF NEBULIZER.
 - CONNECTOR OF NEBULIZER.
 - ADULT MASK LD-N041.
 - CHILD MASK LD-N040.
 - INHALATION TUBE LD-N051.
 - INHALATION MOUTHPIECE.
 - ADULT NASAL NOZZLE LD-N058.
 - CHILD NASAL NOZZLE LD-N059.
 - BAFFLE «A» (YELLOW) LD-N001.
 - BAFFLE "B" (BLUE) LD-N002.
 - BAFFLE «C» (RED) LD-N003.
 - FILTER INHALED.
 - BAG.
 - CONNECTOR OF NEBULIZER.

NEBULIZER THERAPY - WHAT IS IT?
 Nebulizer is a device for formation and spraying of aerosol. The word "nebulizer" is derived from the Latin word "nebula" (fog, cloud) and was first used in 1874 for a device that turns a liquid substance into an aerosol for medical purposes. One of the first portable "aerosol apparatuses" was created by J. Sales-Girons in Paris in 1859. The first nebulizers were used as steam jet energy sources and were applied for inhalation the vapors of resins and antiseptics by tuberculosis patients. Presently, the term "inhaler" is often used instead of "nebulizer". The purpose of the nebulizer therapy is to quickly deliver to the respiratory passages a therapeutic dose of a preparation in aerosol form. Continuous supply of aerosol allows, within several minutes, creating high concentration of a medicine in the upper and lower respiratory passages and lungs, with low probability of any by-effects. Respectively, effective bronchodilation (bronchi expansion) is reached, and the need for hospitalization is eliminated or the hospital stay is reduced. Little Doctor International (S) Pte. Ltd. offers you to use inhaler LD-220C, whose distinctive features are the possibility to use a wide range of medicines, low inhalation solution residual volume, and reliable and simple use. We thank you for your choice.

GENERAL INFORMATION
 Compressor nebulizer LD is designed for treating the diseases of respiratory passages and lungs by medicine solution aerosols. This Instruction Manual is designed to assist the user with safe and effective operation of the Compressor Nebulizer LD. Use this device according to the rules described in this Manual. Operate the Device only as intended. Do not use the Device for any other purposes. Read and understand the whole Instruction Manual. Functionally, the device consists of an air compressor and nebulizer (aerosol formation chamber). The air compressor, on/off power switch and air filter are united in one casing. From the air compressor, the compressed air is directed through a pipe to the nebulizer, where aerosol is formed. For cooling the compressor, air is force-feed into the device's casing.

SAFETY INFORMATION
 To assure the correct use of the product, basic safety measures should always be followed including the warnings and cautions listed in this instruction manual.

WARNING

- For regime of medication shall follow the instructions of your physician or licensed healthcare practitioner.
- Do not cover the compressor with a blanket, towel, or any other type of cover during use. This could result in the compressor overheating or malfunctioning.
- Do not use the device where the device may be exposed to flammable gas or vapors.
- Do not use running/municipal water for inhalations. Use only saline solution.
- Always dispose of any remaining medication in the medication tank after each use. Use fresh medication each time you use the device.
- Do not leave the device or its parts where it will be exposed to extreme temperatures or changes in humidity, such as leaving the device in a vehicle during warm or hot months, or where it will be exposed to direct sunlight.

CAUTION

- Limit the use of the device to 20 minutes at a time, and wait 40 minutes before using the device again.
- Provide close supervision when this device is used by, on, or near infants, children or compromised individuals.
- Do not insert any object into the compressor.
- Make sure that the air filter is clean. If the air filter has changed color or has not been used for 60 days, replace the filter.
- Make sure that the nebulizer kit is correctly assembled, the air filter is properly installed, and the air tube is correctly connected to the compressor and the nebulizer kit. Air may leak from the air tube during use if not securely connected.
- Inspect the compressor (main unit) and the nebulizer parts each time before using the device. Make sure no parts are damaged, the nozzle and air tube are not blocked and the compressor operates normally.
- Do not use the device if the air tube is bent.
- Do not block the air filter cover.
- Do not alter the baffle, the nozzle in the medication tank or any part of a nebulizer kit.
- Do not add more than 10ml of medication to the medication tank.
- Do not operate the device at temperatures greater than 40°C.
- Do not tilt the nebulizer kit so the angle of the kit is greater than 45°.
- Do not shake the nebulizer kit while using the device.
- Do not subject the compressor, or any of the components to strong shocks, such as dropping on the floor.
- This device is approved for human use only.
- Do not disassemble or attempt to repair the device or components.
- Use the device only for its intended use as described in the instruction manual. Do not use attachments not recommended by the manufacturer.
- Dispose of the device, components and optional accessories according to applicable local regulations. Unlawful disposal may cause environmental pollution.
- Make sure that the air tube is securely attached to the compressor (main unit) and nebulizing parts, and does not come loose. Twist the air tube slightly when inserting it into the connectors to avoid the tube disconnecting during use.

RISK OF ELECTRICAL SHOCK

- Do not use the compressor (main unit) and the power cord while they are wet.
- Do not plug or unplug the power cord into the electrical outlet with wet hands.
- Do not immerse the compressor (main unit) in water or other liquid.
- Do not spill water or other liquids on the compressor. These parts are not waterproof. If liquid spills on these parts, please unplug the power cord and wipe off the liquid with gauze or other soft absorbent material immediately.
- Do not use or store the device in humid locations or outdoors. Use the device within the operating temperature and humidity.
- Do not overload power outlets. Plug the device into the appropriate voltage outlet.
- Do not use extension cords. Plug the power cord directly into the electrical outlet.
- Unplug the power cord from the electrical outlet after using the device. Never leave this product unattended when plugged in.
- Unplug the power cord from the electrical outlet before cleaning the device.
- Completely read all of the instructions including the optional accessories before using them.
- Not to touch the ME EQUIPMENT so that it is difficult to operate the disconnection device.
- The power switch is used to isolate the device from the supply mains.
- The direction of movement of the actuator of the supply mains switch is comply with IEC 60447.

MAINTENANCE AND STORAGE

- Keep the device out of the reach of unsupervised infants and children. The device may contain small parts that can be swallowed.
- Do not leave the cleaning solution in the nebulizer parts. Rinse the nebulizer parts with clean hot tap water after disinfecting.
- Wash the nebulizer parts after each use. Dry the parts immediately after washing.
- Do not store the air tube with moisture or medication remaining in the air tube. This could result in infection as a result of bacteria.
- Store the device and the components in a clean, safe location.
- Do not carry or leave the nebulizer with medication in the medication tank.
- Do not place or attempt to dry the device, components or any of the nebulizer parts in a microwave oven.
- Do not wrap the power cord around the compressor (main unit).
- The followings are maintenance and repair which can be taken by operator, or which must be operated by manufacturer or distributor.

SERVICE AND MAINTENANCE	RESPONSIBLE
Change the inhalation tube	Operator
Change the applied part	Operator
Change the air filter	Operator
Clean the surface of the device	Operator
Disinfecting and disinfecting	Operator
All components which need to be repaired or changed by separating the device	Authorized service center of distributor or manufacturer

WARNING:

- Do not modify this equipment without authorization of the manufacturer
- Do not service or maintenance the device while in use with the patient.

CARE AND DISPOSAL

- Regularly clean the device and all accessories. It is recommended that all the accessories should be wiped with a 3% solution of hydrogen peroxide with addition of 0.5% solution of a detergent (for example, a laundry powder). After that, the nebulizer should be washed by a rich jet of water. The mouthpieces and nose nozzles may be treated by boiling for 10 minutes or autoclaving at 150°C. After the treatment, wipe dry all parts of the device with a soft cloth.
- Regularly check the filter for dirt and replace it when needed. FILTER REPLACEMENT IS RECOMMENDED AT

LEAST ONE YEAR.

- Repair the Device only in authorized organizations.
- On expiration of the warranted service life apply from time to time to authorized repair organizations to check the technical condition of the Device. Dispose of time to authorized repair organizations to check the technical condition of the Device. Dispose of the Device and its components according to the application local regulations. No special requirements to disposal of this Device are defined by the manufacturer.

USING THE DEVICE PREPARING FOR INHALATION
 IMPORTANT: Before using the appliance for the first time it is necessary to make a full cleaning, as described in last paragraph «SAFETY INFORMATION».

- Place the nebulizer in front of you on the table. Make sure the device is turned off (power switch is in position «0»), and the power cord is not plugged into the mains.
- Remove the top of the nebulizer by turning it counterclockwise (Fig.2.a).
- Set the desired baffle.
- Factory installed baffle inside the nebulizer is baffle «B» (Blue), which is effective to affect the entire respiratory tract. For a more effective impact on the upper respiratory tract, set instead of the blue baffle, red baffle «C». For a more effective impact on the lower respiratory tract - baffle «A» yellow color, which consists of two parts (Fig.2.b, figures indicate the order of assembly). Graphics of the differential particle size distribution by mass for different baffles are shown in Fig. 3. For effective delivery medicines to a specific site of respiratory tract, need to use the appropriate baffle.
- Fill the bottom of the nebulizer inhalation solution. The dosage should not exceed the recommended by your doctor. The number of nebulizer solution is determined by the scratches on the case. The maximum reservoir volume of 10 ml.
- Attach the nebuliser at the top, turning it clockwise until it stops.
- Depending on the type of inhalation, using either a mouthpiece or nozzle or mask.
- Connect one end of the inhalation tube to the fitting of compressor and others - to the fitting of nebulizer.
- Turn on the nebulizer, switching the power switch in position «0». NEBULIZER IS READY FOR INHALATION.

PERFORMING THE INHALATION
 The length of one treatment session should not exceed 20 minutes. Consult your attending physician about the length of the inhalation procedure. You should always be calm and relaxed during the inhalation. Breathing should be slow and deep, so that the preparation could fill the lungs well and reach the deep portions of bronchi. Briefly hold your breath, and then exhale slowly. Do not attempt to breathe too rapidly. Make pauses if you feel that you need it.

BREATH-ACTUATED NEBULIZER (fig.4)
 The special design of the nebulizer in the form of chambers connected in a certain manner provides different ways of air streams during inhaling and exhaling. It allows obtaining the air stream with greatest aerosol concentration when inhaling and reducing aerosol loss when exhaling. The effectiveness of inhalation using the breath-actuated nebulizer is increased significantly.

COMPLETING THE INHALATION
 When the inhalation solution is used up and the inhalation time recommended by the doctor has expired, turn the device off by putting the tumer in «0» position and unplug it. After inhalation, breathe fresh air for some time for better treatment effect. After each application of the device, the residual preparation should be removed out of it. Clean and wash the device as described in last paragraph «SAFETY INFORMATION».

TECHNICAL SPECIFICATIONS

Model	LD-220C	Max. operating time	20 minutes
Type	Compressor Nebulizer	Operation mode	20 minutes on, 40 minutes off
Rating	AC 230V, 50Hz	Operating temperature, humidity/atmospheric pressure	10°C – 35°C 15%– 80% 86 – 106 kPa
Power consumption	190VA	Transport and storage temperature and humidity/atmospheric pressure	-20°C – 40°C 15%– 95% 50 – 106 kPa
Extreme Pressure	210kPa–400kPa	Net weight (compressor)	~1215 g
Free Flow Range	≥7L/min	Size (compressor)	147 x 95 x 137 mm
Nebulizing Pressure	60kPa–180kPa	Net Weight (nebulizer)	~1215 g
Noise level, approx.	~55 dB*	Pollution Degree	Degrees 2
Max capacity of nebulizer kit	10ml	Overvoltage Category	Category II
Nebulization rate, approximately	Inhalation baffle «A» LD-N001 0.3 ml/min.* Inhalation baffle «B» LD-N002 0.4 ml/min.* Inhalation baffle «C» LD-N003 0.5 ml/min.*	Year and month the manufacture is given in the bottom of the unit body in a serial number after symbol "A"	
Particle size (MMAD)	Inhalation baffle «A» LD-N001 3.5 μm.* Inhalation baffle «B» LD-N002 4.0 μm.* Inhalation baffle «C» LD-N003 5.0 μm.*		

SYMBOL EXPLANATION:

- CE Compliance with the Directive MDR (EU) 2017/745
- Important: Read the instruction
- Authorized Representative in the EU
- Manufacturer
- Keep from moisture
- Protection class II
- BF type product
- Storage/operation temperature, humidity and atmospheric pressure limitation
- Serial number
- Medical device
- Unique Device Identification
- Importer

Revision date of the present Manual is indicated on the last page as XX-XXXX-YYMM-NN, where YY is the year, MM is the month and NN is the number of revision. Technical characteristics may be changed without preliminary notification to improve the operation and quality of the product.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

TROUBLESHOOTING TIPS
 No power on device when the power switch is on:

- Turn the power switch off. Plug the power plug into an electrical outlet. Turn the device on.

No nebulization or low nebulization rate when the power is on:

- Add the correct amount of prescribed medication to the medication cup.
- Make sure the nebulizer kit is correctly assembled and the inhalation accessory is correctly attached.
- Hold the nebulizer kit correctly. Do not tilt the nebulizer kit so the angle of the kit is greater than 45 degrees.
- Make sure the air tube is correctly attached to the compressor and the nebulizer kit.
- Make sure the air tube is not folded, kinked or bent. Inspect the air tube for any damage. Replace the air tube if damaged. The device is very hot:
- Do not cover the compressor with any type of cover during use. Turn the device off. Wait 40 minutes before using the device again.

WARRANTY
 The purchased device is warranted for 3 years from the date of sale. Faults detected during the warranty period will be repaired free of charge by a specialized service center within 21 business days from the date of device delivery to the service center. The costs for the defective device transportation are borne by Little Doctor Europe Sp. z o.o. (door-to-door service), except for faults that are not covered by the warranty. The party submitting the claim must properly prepare the device for shipment, protecting it from damage during transportation and enclose a completed warranty card (with the seller's stamp, date of sale and seller's signature) and proof of purchase (invoice or receipt). The warranty does not cover damages resulting from inadequate protection of the device during transportation. The warranty period is extended for the time the device is under repair. Post-warranty repair of the device is performed at the user's expense. The warranty does not cover the following: consumables and accessories for individual use (tubes, nebulizers, masks, mouthpieces, etc.); products with mechanical, thermal or chemical damage, or damage caused by violations of the rules specified in the operating instructions; products with signs of modification, opening and/or repair by an unauthorized service center (by an individual); products that have been connected to unrecommended power supplies or used with unrecommended consumables; products that have been operated with non-remedied defects; products that are defective through no fault of the manufacturer, such as: actions of third parties, natural phenomena and natural disasters, animals, insects, foreign objects or liquids entering the device, etc.; internal and external contamination, scratches, cracks, scrapes and other mechanical damage to the equipment occurring during its operation; products with illegible or altered serial numbers or a damaged warranty label; batteries, glass parts, light bulbs, as well as operations such as adjusting, cleaning and other maintenance of the product. The buyer has the right to replace the device with a serviceable one if an authorized service center finds a manufacturing defect that cannot be eliminated or after 3 unsuccessful repairs of the same malfunction.

MANUFACTURER AND AUTHORIZED REPRESENTATIVES
 This device manufacturing is certified according to international standard ISO 13485. Device comply with the requirements of European Directive MDR (EU) 2017/745, international standards EN 1041, EN 10601-1, EN 10601-1-2, ISO 14971.

Manufactured under control and for Little Doctor International (S) Pte. Ltd. (7500A BEACH ROAD #11-313 THE PLAZA SINGAPORE 199591. Postal Address: Senang Central PO Box 960 Singapore 915411).
Manufacturer: Little Doctor Electronic (Nantong) Co. Ltd. (No. 8, Tongyang Road Economic & Technical Development Area, 226010 Nantong, Jiangsu, People's Republic of China).
 Little Doctor Europe Sp. z o.o., 576 Zawiła Street Krakow 30-390 Poland.
 Little Doctor Europe Sp. z o.o., 576 Zawiła Street Krakow 30-390 Poland.

IMPORTANT INFORMATION REGARDING ELECTRO MAGNETIC COMPATIBILITY (EMC)
 This Device manufactured by Little Doctor Electronic (Nantong) Co. Ltd. conforms to IEC-C60601-1-2:2014 Electro Magnetic Compatibility (EMC) standard. Further documentation in accordance with this EMC Standard is available at Representative in EU Little Doctor Europe Sp. z o.o. at the address mentioned in this instruction manual or at www.littledoctor.sg. Refer to the EMC information for Device on the website.

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