Meriden Comfort SOFT

Electric Heating pad

Elektrisches Heizkissen

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CONGRATULATIONS!

Thank you for your purchase of Meriden Comfort SOFT Electric Heating pad. Please read this manual carefully and keep it for future reference.

IMPORTANT SAFETY INSTRUCTIONS. READ CAREFULLY AND KEEP FOR FUTURE REFERENCE.

Do not use folded or rucked.

Do not insert pins.

1. This product can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they are being supervised or have been instructed concerning the usage of the product in a safe way and if they understand the hazards involved.
2. Children shall not play with the product.
3. Cleaning and user maintenance shall not be carried out by children unless they are older than 8 and supervised.
4. This product can be used by children older than 3 years and younger than 8 years under supervision and with the control always set to minimum temperature value.
5. If the supply cord is damaged, it must be replaced by the manufacturer, its service agent or a similarly qualified person in order to avoid hazards.
6. When the product is likely to be used for an extended period, controls should be adjusted to a setting recommended for continuous use.
7. The product is only to be used with the control units types (MCS-301) that are marked on the product.
8. When storing the product, allow it to cool down before folding.
9. Do not crease the product by placing items on top of it during storage.
10. Examine the product frequently for signs of wear or damage. If there are such signs or if the product has been misused, return it to the supplier prior to any further use.
11. This product is not intended for medical use in hospitals.
12. Do not use the product when wet.
13. The product must not to be used by persons insensitive to heat and other very vulnerable persons who are unable to react to overheating.

14. Children under the age of three are not to use this appliance due to their inability to react to overheating.
15. The cable and control unit of the product can give rise to risks of entanglement, strangulation, tripping or treading if not correctly arranged. The user shall make sure that excess ties and electric cords shall be arranged in a safe way.

16. Read this manual thoroughly before first use thoroughly. It contains important information for your safety as well as for the use and maintenance of the equipment.
17. Keep this manual for future reference and pass it on with the product.
18. Pay particular attention to the warnings on the product and in this manual.
19. Use the product only for its intended purpose. Improper use may lead to hazards.
20. In the event of improper use or wrong operation, we do not accept any liability for resulting damages.
21. The use of accessories and components not expressly recommended by the manufacturer may result in injuries or damages and voids the warranty.
22. Never leave the product unsupervised while it is switched on.
23. Do not place any objects on top of the product.
24. Never immerse electrical parts of the product in water or any other fluid during cleaning or operation. Never hold the product under running water.
25. Do not under any circumstances attempt to repair the product yourself. For servicing and repairs, please consult an authorised service agent.
26. The mains voltage must match the information on the rating label of the product.
27. The power plug is used as the disconnect device. Please be sure that the power plug is in operable condition.
28. Only unplugging from the power socket can separate the product completely from the mains.
29. Avoid damages to the power cord that may be caused by kinks or contact with sharp corners.
30. Keep this product including power cord and power plug away from all heat sources such as ovens, hot plates and other heat-producing products/objects.
31. Only unplug from the power socket by the plug itself. Do not pull on the cord.

32. Unplug the product ...
 - when the product is not being used,
 - in case of malfunctions,
 - before attaching / detaching accessories and
 - before cleaning the product.
33. If the power cord or the casing of the product are damaged or if the product has been dropped, the product may not be used before it is examined by an authorised service agent.
34. Use this product in dry indoor areas only and never outdoors.
35. Protect this product from extreme heat and cold as well as from direct sunlight and dust.
36. In order to reduce the risk of fire and electric shock, do not expose the product to fire (fireplace, grill, candles, cigarettes and so on) or water (water drops, water splashes, vases, tubs and so on).
37. Warning! Do not use this product near water or heavy moisture; e.g. in a damp basement or beside a swimming pool or bathtub.
38. Use this product in temperate latitudes only. Avoid tropical and especially humid climatic regions.
39. Protect this product from all moisture including spray water and humidity.
40. This product is not suitable for commercial use. It is designed for household use only.
41. The product is only to be used with the types that are marked on the product.
42. The switch or control unit must not be allowed to get wet during washing and during drying the cord must be positioned in a way, that ensures, that water does not flow into the switch or control.
43. Do not fall asleep while the product is switched on.
44. Prolonged use of the product at a high setting could lead to skin burns.
45. Follow the instructions in the chapter Cleaning and care.

INTENDED USE

This product is designed for warming the human body only. Any other use may result in damage to product or injuries. The Imtron GmbH assumes no liability for damage to the product, for property damage or for personal injury due to careless or improper usage of the product, or usage of product which does not meet the manufacturer's specified purpose.

CONTROLS AND COMPONENTS

- A. Heating pad
- B. Connector jack
- C. Power cord with plug
- D. Operation light
- E. Switch (0 1 2 3)

DELIVERY

- 1 x Heating pad
- 1 x Instruction manual
- 1 x Warranty card

BEFORE FIRST TIME USE

- Remove the product and accessories carefully from the original packaging. It is recommended to keep the original packaging for storage. If you wish to dispose of the original packaging, please observe applicable legal provisions. Should you have any questions regarding proper disposal, contact your local waste management center.
- Inspect the delivery contents for completeness and damages. Should the delivery contents be incomplete or damaged, contact your sales outlet immediately.
- Carefully check the heating pad before each use for wear or damage.
- Do not use if you notice wear, damage or signs of improper use on the heating pad
- the switch or the cables or if the product is not working.

USING THE HEATING PAD

Fig. 1 Connect the pad connector with the pad.

Fig. 2 Connect the power plug to a suitable power outlet.

Fig. 3 Switch on the heating pad.

Note

Please note there is a feature to show whether the pad connector is well connected to the pad. When the display keeps flashing, please firstly switch off the appliance and pull out the power plug from the socket outlet, then check whether the connector is well connected to the pad, or reconnect the connector to the pad, then plug in and switch on the product again.

Fig. 4 Set the desired temperature.

Setting	Function
0	Off
1	Minimum heat
2	Medium heat
3	Maximum heat

Note:

- If the heating pad is too warm, please reduce the temperature setting.
- Setting 1 is recommended for continuous use. This heating pad is not intended for use in hospitals or for commercial applications.
- It may not be used on infants, small children, persons insensitive to heat, helpless persons or animals.
- Do not fall asleep while the heat pad is switched on.
- Prolonged use of the heating pad at a high setting could lead to skin burns.

AUTOMATIC SWITCH-OFF FUNCTION

Fig. 5 The heating pad switches off automatically after approx. 90 minutes continuous operation. Switch off the heating pad (0) or unplug the power plug from the power outlet if you are not going to use the heating pad after the automatic switch-off.

CLEANING AND CARE

- Warning! Before cleaning, always let the product cool down and disconnect it from the mains.
- Warning! When cleaning, never use solvents or abrasive materials, hard brushes, metallic or sharp objects.
- Warning! When cleaning, never use solvents or abrasive materials, hard brushes, metallic or sharp objects. Solvents are harmful to human health and can attack plastic parts, while abrasive cleaning mechanisms and tools may scratch the surface(s)

Fig. 6 Unplug the power plug from the power outlet. Detach the connector from the heating pad. Remove smaller stains with soft brush.
Clean the textile cover according to the laundry care symbols on the label.

	Machine washable up to 30° C		Do not dry clean.
	Do not tumble dry.		Do not bleach.
	Do not iron.		

Note

- Do not spin or machine dry, chemically clean, wring out, press or iron the heating pad.
- Only use the heating pad again when it is completely dry.

SPECIFICATION

Rated voltage: 220 - 240 V~, 50 Hz

Rated power: 100 W

Protection class: II

DISPOSAL

Do not dispose of this product as unsorted municipal waste. Return it to a designated collection point for the recycling of WEEE. By doing so, you will help to conserve resources and protect the environment. Contact your retailer or local authorities for more information.

Meriden Sp. z o.o.
ul. Piastów 96
42-290 Blachownia, Poland
meriden@meriden.pl
www.meridenglobal.com

www.facebook.com/meridenglobal
www.instagram.com/meridenglobal
Designed by Meriden in Poland
Made in P.R.C.

www.meriden.pl