Meriden Sp. z o.o.
ul. Piastów 96,
42-290 Blachownia, Polska
NIP: 9492217625
meriden@meriden.pl
www.meridenglobal.com

www.facebook.com/meridenglobal
www.instagram.com/meridenglobal

Designed by Meriden in Poland
Made in P.R.C.

Model: PG-IRT1603
Date: 2022-09-09

Infrared Ear/Forehead Thermometer

Thanks for buying and using this product, please read this manual carefully before use.

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Shenzhen Pango Medical Electronics Co., Ltd

Main Site: No.25 1st Industry Zone, Fenghuang Road, Xikeng Village, Henggang Town, Longgang District, Shenzhen, Guangdong China

Additional site1: 2-4 Floor, No.5 Shanzhuang Rd., Xikeng Village, Henggang Town, Longgang District, Shenzhen City, Guangdong Province, China
P.C.: 518155

Web: <http://www.pan-go.com>

Authorised representative: Lotus NL B.V.

Address: Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherland

Tel: +31645171879 (English), +31626669008 (Dutch)

REP

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PG-IRT1603

Catalogue

Safety signs.....	2
Product introduction	3
User's safety	4
Temperature taking hints.....	5
Product description	6
The screen display marks description	7
Battery installation procedure.....	7
Basic functions guideline.....	8
Temperature taking methods.....	8
Care and cleaning	10
Frequently Asked Questions and Solutions	11
Troubleshooting	12
Product specification	13

EN

Safety signs

- ◆ The warning signs and illustrations indicated in the manual are intended to enable you to use the product safely and correctly and to prevent any harm to you and others
- ◆ Warning signs, illustrations and their meanings are as follows.

 **Caution: please refer to the below table**

The key	
	Points to observe.
	General warning.
	General prohibition.
	No disassembly allowed.
	Equipment with type BF applied parts.
	Water and moisture resistance.
	Refer to the manual before use.
	Packaging of the unit may comply with the requirements of environment protection.
	Material and the product itself is made of renewable material allowing recycling of the product to the benefit of the environment and the planet.
	Do not dispose of the product in the household waste.
IP Classification: IP22	

EN

Warning

<ul style="list-style-type: none"> ● Measured values are to be consulted only with a qualified practitioner. ● The product is to be used only for measuring of human body temperature. ● Please do not use this product for purposes other than human body temperature measurement. ● No mobile phones are allowed to be used in close proximity of the product. ● Please do not use equipment that generates electromagnetic fields in close proximity of the product. 	
<ul style="list-style-type: none"> ● Do not disassemble or repair the device using spare parts on your own. ● Do not attempt to bend or stretch the device ● Protect the product from dropping or a severe impact. 	

Warning

- The unit is intended for measuring of a human body temperature only and cannot be used in disease diagnosis, emergency or for continuous measurement in a surgery.
- Keep the product out of reach of children. Children up to 12 years old should use the product under supervision of their parents.
- Patients cannot attempt disease diagnosis or treat themselves on the basis of temperature measurement result. It is strongly advised that a professional consult is sought.
- Children under the age of 12 are prohibited from using the product.
- The product should not be used by people suffering from external otitis, tympanitis or other ear diseases.

If the product is used or stored beyond specified temperature and humidity range, it may fail to reach original performance provided in the specification.
Operational environment: temperature: from +10°C~+40°C, humidity: from 15%RH~93%RH
Storage environment: temperature: from -25°C ~+55°C, humidity: from 0%RH~93%RH

Product introduction

Intended use: Infrared Forehead Thermometer is intended to measure human body temperature by measuring in ear canal or forehead.

Scope of application: It is suitable for displaying body temperature of the measured object by measuring the heat radiation in ear canal or forehead.

Features:

1. Non-touch infrared measurement of ear/forehead temperature.
2. Multiple colors and backlight display: white, green, orange and red.
3. Possibility to store up to 9 temperature measurements.
4. The switch between degree Fahrenheit °F and degree Celsius °C. (original setting: degree Celsius °C)
5. Instant measurement within 1 second.
6. Convenient and economical design without the use of an earmuff cuts on subsequent use costs.
7. Sound on/off function.
8. If idle for more than 30 seconds, it turns off automatically.

Hint:

Measurement result are only a reference point and cannot, at any circumstances, replace a medical diagnosis of a qualified practitioner. Should you have any questions regarding the measurement result, please seek a doctor's consult.

User's safety

Warning

1. Always seek a professional consult as it is dangerous to diagnose and treat on the basis of measurement result.
 - ◆ self diagnosis may lead to worsening of a patient condition.
2. Please do not attempt to touch or blow on the infrared sensor.
 - ◆ Dirty or damaged infrared sensor may cause abnormal measurement results.
3. If there is a difference in temperature between the storage site and measurement site, please place the thermometer in a room temperature (at measurement site) for about 30 minutes before a measurement.
 - ◆ Difference in temperatures may cause incorrect measurement results.
4. Keep the product out of reach of children.
 - ◆ Children may hurt themselves while handling or attempting to measure temperature. If a child swallows the battery or the transparent cover, contact a doctor immediately.
5. Do not attempt to measure body temperature near the air conditioning unit.
 - ◆ Avoid affecting measurement accuracy.
6. Before and after every use, wipe the probe surface with a cotton swab dipped in 95% absolute alcohol (if you notice stains, fog or water on the infrared sensor, wipe it gently with a cotton swab dipped in 95% anhydrous alcohol).
 - ◆ Cleaning of infrared sensor with the use of toilet paper or facial tissue may leave scratches which in turn may affect measurement result.
 - ◆ Avoid cross-infection and affecting accuracy of the measurement.
7. Protect the product from suffering mechanical damage.
 - ◆ Damaged product may cause inaccurate measurement results.
8. If the product gets accidentally immersed in water or splashed, be sure to dry the product completely. Water on the sensor surface should be removed with a clean cotton swab.
 - ◆ Avoid causing accidents and affecting measurement accuracy.

Caution:

1. People suffering from external otitis, tympanitis and other ear diseases should not attempt to use the product.
 - ◆ It is possible for the affected area to worsen.
2. The product should not be used directly after swimming, bathing or when one's ears are wet.
 - ◆ Measured temperature might be low.
3. Do not expose the battery to risk of damage.
 - ◆ The battery might suffer damage and break.
4. When measuring temperature in ear canal the product must operate in ear temperature mode.
 - ◆ Measured value might be incorrect.

Suggestions

1. When referring measured values to a doctor, please highlight that the temperature was taken with the use of an ear thermometer.
2. Do not deliberately expose the product to risks of impacting, falling, trampling or shaking.
3. Do not attempt to disassemble, repair or modify the product.
4. Do not allow liquids (alcohol, water drops, hot water etc.) to enter the product body as it is not water resistant.
5. The product must be kept clean and stored in a dry place.
6. If any problems with functioning of the product are noticed, please contact sales department. Do not attempt to repair the product on your own.
7. Do not use the product in environment where electromagnetic interferences occur.
8. At the end of product useful life, please dispose of it according to local laws and regulations.

Temperature taking hints

Comparison of different measurement methods.

Measured values are different if different measurement methods are used. The WHO provides normal human body temperature reference values, please see below table about the specific temperature difference.

Measurement method	Normal body temperature
Rectal temperature	36.6°C-38°C

Oral temperature	35.5°C~37.5°C
Axillary temperature	34.7°C~37.3°C
Cochlear temperature	35.8°C~38°C
Oral temperature	35.5°C~37.5°C (PG-IRT1603 measured value)

Changes in human body temperature

Human belongs to species with a constant temperature. Body temperature is basically constant, however it is not totally changeless. Human body temperature changes during a day and can be observed as follows:

At night

Lowest Body temperature is lowest due to sleep and decreased activity.

In the morning

Higher Moving from warm bed to lower temperature of the room in the morning causes body muscles to contract and produce heat.

At noon

Highest human body reaches the highest temperature after lunch causing body to adjust naturally.


Three or four o'clock in the afternoon

Lower Due to physical exertion, decreased sugar levels in blood.


In the evening

Lowest Due to the sun going down, decreased room temperature.

Product description

The screen display marks description


Battery installation procedure

Installing the battery:


1. Press the battery cover to remove it, it will pop out automatically.
2. Prepare two AAA batteries of 1.5V (it is recommended to use alkaline batteries). Insert the batteries with the correct polarity into the battery compartment as indicated (+/-).

Press the battery cover to remove it

Low battery caution:

When the steady battery icon shows on the display together with „LO“, it means the device cannot work properly and the batteries must be replaced immediately.


Hint


- Remove the battery from the product if it is not required for extended periods of time in order to avoid damage to the thermometer resulting from a leaking battery.
- It is recommended to use alkaline batteries 
- To protect the environment, dispose of empty batteries at appropriate collection sites according to national or local regulations

Basic functions guideline

Sound on/off mode:


1. In the power-on mode, press „mem“ button to turn the sound on or off.
2. After pressing the „mem“ button, the icon  will appear, together with a short sound signal, which indicates that the sound is in ON mode.
3. After pressing the „mem“ button  again, the icon  will appear, which indicates that the sound is in silent mode.

Changing the temperature scale from °F to °C

In the power-off mode, press „mem“ button for six seconds- to set preferred temperature scale: Fahrenheit or Celsius (°F or °C). The device will switch off automatically after 8 seconds or it can be switched off immediately by pressing the button .

Memory function

In the power-off mode, press „mem“ button- the device has access to the last nine measurement results in order from the latest to the oldest (see the picture below). The device will switch off automatically after 30 seconds of no use, or it can be switched off immediately by pressing the button .

Backlight function

When the measurement result indicates the temperature <34°C, the device will show LO sign with red backlight.

When the measurement result indicates the temperature it means normal body temperature and is shown with green backlight.


When the measurement result indicates the temperature it means slightly raised body temperature and is shown with orange backlight.

When the measurement result indicates the temperature, it means high body temperature and is shown with red backlight.


When the measurement result indicates the temperature >43.0 °C, the device will show HI sign with red backlight.

Hint: This function is for reference only.


Temperature taking methods

Ear temperature

Ear temperature


Forehead temperature

1. Ear temperature measurement

Hint:


For the ear temperature measurement, remove the head cover intended for forehead measurement.

1.1 After removing head cover, press button **ⓘ** to switch on the device.

When the ear temperature measurement mode is on, the last measurement result will be displayed. Put the thermometer in the ear and press the button **ⓘ** to measure the temperature.

1.2 After 1 second, the measurement result will be displayed.

Flashing LED screen indicates the measurement preparation


Measurement result

Note: Until the sound signal is heard, the measurement is not completed. Only after the beep sound appears, can the probe be taken out from the ear. (If the sound indication has been turned off, the sound signal will not appear.)

2. Forehead temperature measurement


Hint:

Install the head cover intended for forehead temperature measurement.

After installing the head cover, press button **ⓘ** to switch on the device. When the forehead temperature measurement mode is on, the last measurement result will be displayed. Aim the thermometer at the central part of the forehead and keep it within the distance of 0-5 mm from the forehead. Press the button **ⓘ** to measure the temperature.

After 1 second, the measurement result will be displayed.

Flashing LED screen indicates the measurement preparation

Measurement result


Note: Until the sound signal is heard, the measurement is not completed. Only after the beep sound appears, can the probe be taken out from the ear. (If the sound indication has been turned off, the sound signal will not appear.)

3. Object temperature measurement

Hint: The installation of head cover used for forehead temperature measurement has no influence on the result of object temperature measurement.

3.1 Press the button **ⓘ** for 6 seconds to enter the object temperature measurement mode. Aim the thermometer at the object and press the button **ⓘ** to measure the temperature.

3.2 After 1 second, the measurement result will be displayed.



Flashing LED screen indicates the measurement preparation

Measurement result

Note: Until the sound signal is heard, the measurement is not completed. Only after the beep sound appears, can the probe be taken out from the ear. (If the sound indication has been turned off, the sound signal will not appear.)

Care and cleaning

1. Use cotton swab moistened with absolute alcohol to clean the measuring probe after each use. It will help to prevent cross-infection. (Note: Never immerse the device in water.)
2. Use a soft and dry cloth to wipe the thermometer casing. Do not use any agents which could cause scratches on the device. Never clean the device with water or other cleaning liquids.

Hint:

Why cleaning after each measurement is needed?

Infrared thermometer uses highly sensitive technique to measure temperature of the target object. Not only will ear wax or dust affect the measurement accuracy but also make bacterial infection more likely. It is recommended to clean the product as shown in the above illustration after every use.

Frequently asked questions and solutions

Displayed prompt	Reason	Solution
Hi	When the target object temperature is higher than the measurement range of ear temperature mode set to maximum of 43.0°C, LCD screen will indicate "HI" prompt.	1. (When the measurement probe is not properly placed in the ear canal during the measurement process or the distance is too far, the measurement result might be low.)
Lo	When the target object temperature is lower than the measurement range of ear temperature mode set to minimum of 34.0°C, LCD screen will indicate "Lo" prompt.	2. When measurement probe is dirty, the measurement value might be low. Use a cotton swab dipped in alcohol to wipe the measurement probe.
Er.H	The product has an upper limit of operating temperature of 40°C. When the ambient temperature exceeds this temperature point, the LCD screen will display an error message „Er.H“.	The environment temperature cannot be greater than 40°C while operating the product.
Er.L	This product has a lower limit of operating temperature of 10°C. When the ambient temperature exceeds this temperature point, LCD screen will display an error message „Er.L“.	The environment temperature cannot be lower than 10°C while operating the product.
Err	When the environment temperature changes rapidly for about 5 °C during measurement in object temperature mode, an „Err “ prompt will be indicated then the product shutdowns automatically.	When message „Err“ is displayed, leave the product in current environment for 30 minutes before next measurement.

Should any operation go wrong, If the sound is on, a short beep will be heard.

Troubleshooting

Event	Reason	Solution
The display does not turn on when the power is switched on.	Battery is exhausted,	Replace the battery
Measurement temperature is low.	Battery polarity is wrong.	Make sure that the battery polarity is the same as in the battery case.
	Measurement position is not correct.	Follow the instructions while measuring the temperature.
	Ear canal or the sensor itself is dirty.	Remove the dirt before measurement.
Substantial temperature fluctuations during continuous measurement.	Measurement intervals are too short.	The interval for each measurement should be above 10 seconds.

Product specification

Product name: Infrared Ear/Forehead Thermometer

Model number: PG-IRT1603

Product dimensions: 31×175×72mm

Product weight: about 77g (excluding battery)

Measuring range: 34.0°C-43.0°C.(93.2-109.4°F)

Object temperature: 0- 93.2°C (32-199.7°F)

Resolution ratio: 0.1°C/°F

Measurement Location: laboratory

Accuracy: (35.0°C~42.0°C)±0.2°C (95.0~107.6°F)
±0.4°F other temperature±0.3°C.

Operational temperature: 10.0°C~40.0°C(50.0°F~104.0°F)
relative maximum humidity 15%RH~93%RH

Atmospheric pressure: 70kPa~106kPa

Transportation/storage temperature: -25~55°C (-13°F~131°F),
relative maximum humidity 0%RH~93%RH

Atmospheric pressure: 50kPa~106kPa

Display screen: LCD display screen, 4 digits and special icons.

Sound: when switched on a short beep will be heard indicating that the product is ready to measure. When measurement is finished a long beep will be heard.

System error or fault: three consecutive short beeps.

Fever alert: ten consecutive short beeps.

Memory: up to nine measurements can be stored.

Automatic shutdown: if no operation for 30 seconds, the thermometer will shutdown automatically.

Battery: two 1.5V AAA batteries (alkaline batteries are recommended).

Period of use: five years

1603 Forehead mode:

Clinical bias, Dcb: 0.078

Limits of Agreement, LA 0.243

Clinical Repeatability, or 0.069

Referenced body site: forehead

Measuring site: forehead

Packing parts list

1. Main body
2. User manual

Appendix 1 Guidance and Manufacturer Declaration Tables

Guidance and manufacturer's declaration –electromagnetic emissions		
Model PG-IRT1603 is intended for use in the electromagnetic environment specified below. The customer or the user of Model PG-IRT1603 Infrared Thermometer should assure that it is used in such an environment.		
Emissions	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	Thermometer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	Thermometer is intended to be used at home and it is powered by DC 3V
Harmonic emissions IEC 61000-3-2	N.A.	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	N.A.	

Guidance and manufacturer's declaration – electromagnetic immunity.

Model PG-IRT1603 Infrared Thermometer is intended for use in electromagnetic environment specified below. The customer or the user of Model PG-IRT1603 Infrared Thermometer should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD)IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15KV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 KV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m, 50/60Hz	30 A/m, 50/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE UT is the AC mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration electromagnetic immunity

The is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-IRT1603 Infrared Thermometer should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms 150 kHz to 80 MHz outside ISM bands	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the Model PG-IRT1603 Infrared Thermometer, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m	<p>Recommended separation distance</p> $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 80\text{MHz to } 2.7\text{GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range b Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Model PG-IRT1603 Infrared Thermometer is used exceeds the applicable RF compliance level above, the Model PG-IRT1603 Infrared Thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Model PG-IRT1603 Infrared Thermometer.

Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Model PG-IRT1603 Infrared Thermometer

The Model PG-IRT1603 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Model PG-IRT1603 Infrared Thermometer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Model PG-IRT1603 Infrared Thermometer as recommended below, according to the maximum output power of the communications equipment. Model PG-IRT1603 Infrared Thermometer

Rated maximum output of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.7 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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